


HIMEDIA LABORATORIES PVT.LTD

INDIA



WEBINAR TOPIC

**HIENCAP™
MEDIA**

**GRANULATED
MEDIA**

**CHROMOGENIC
MEDIA**

HIMEDIA®

For life is precious

HiEncap™ Culture Media



ADVANCEMENT IN CULTURE MEDIA FORMS



19th Century-
Media
preparation with
minced meat
and other
ingredients



20th century-
dehydrated culture
media



21st century-
Encapsulated
culture media



HiEncap™ Culture Media

Features

- User friendly
- Simple boiling required
- Accurate and reliable

Advantages

- Premeasured media
- Saves time ,avoid hassles of weighing and wastage
- Ready to be autoclaved
- pH adjusted
- Prevents generation of aerosols

Applications


- Clinical, Food, Water, Dairy, Industries, Cosmetics, Pharma, Molecular Biology



HiEncap™ Culture Media

Features

- Dehydrated media enclosed in a gelatin capsule free of TSE/BSE risks
- Three easy steps :
Simply suspend in water, autoclave and use
- Premeasured quantity stored in each capsule
(250 ml, 500ml, 1000ml)



HiEncap™ Culture Media Capsules




Range of HiEncap™ Culture Media for:
100 ml, 250 ml, 500 ml & 1000 ml

HiEncap™ Media Preparation




Simple requirement, HiEncap™ Culture
Media & Water for preparation


Step A - Suspending HiEncap™ Culture Media in Distilled Water / Purified Water



Preparing 100 ml Media



Preparing 250 ml Media



Preparing 500 ml Media



HiEncap™ Culture media suspended

Step B - Heating HiEncap™ Culture Media to Dissolve & Digest the Medium (100°C or BWB or Microwave Oven)



Heating 100 ml Media



Heating 250 ml Media



Heating 500 ml Media



Heating HiEncap™ Culture media in
Microwave Oven


Step C - Pre Heating Facilitates Digestion of HiEncap™ Culture Media



Digestion of 100 ml Media



Digestion of 250 ml Media





Digestion of 500 ml Media





Digestion of HiEncap™ Culture Media

Step D - Homogeneous Medium Ready for Autoclaving



Step E - Homogeneous Medium after Autoclaving



HiEncap™ Culture Media

Code	Product
EC073	HiEncap™ Blood Agar Base
EC211	HiEncap™ Brain Heart Infusion Agar
EC210	HiEncap™ Brain Heart Infusion Broth
EC614/EC1494I	HiEncap™ Buffered Peptone Water
EC557	HiEncap™ Luria Agar
EC1726	HiEncap™ Luria Agar Base, Miller's Modification
EC1151	HiEncap™ Luria Bertani Agar, Miller
EC1245	HiEncap™ Luria Bertani Broth, Miller
EC575	HiEncap™ Luria Broth
EC1725	HiEncap™ Luria Broth Base, Miller's Modification
EC081	HiEncap™ MacConkey Agar w/0.15%Bile Salts, CV and NaCL
EC082A	HiEncap™ MacConkey Agar w/o CV, NaCL w/0.5% Bile Salts
EC083	HiEncap™ MacConkey Broth Purple w/ BCP
EC007	HiEncap™ MacConkey Broth w/ Neutral Red
EC173	HiEncap™ Mueller Hinton Agar
EC1084	HiEncap™ Mueller Hinton Agar No. 2
EC391	HiEncap™ Mueller Hinton Broth
EC001	HiEncap™ Nutrient Agar
EC002	HiEncap™ Nutrient Broth
ECG020	HiEncap™ NZYM Growth Medium
EC091	HiEncap™ Plate Count Agar

Code	Product
EC096	HiEncap™ Potato Dextrose Agar
EC403	HiEncap™ Potato Dextrose Broth
EC063	HiEncap™ Sabouraud Dextrose Agar
EC033	HiEncap™ Sabouraud Dextrose Broth
ECG016	HiEncap™ SOB Growth Agar
ECG014	HiEncap™ SOB Growth Medium
EC1252	HiEncap™ SOB Medium
EC1379	HiEncap™ SOC Broth
ECG015	HiEncap™ SOC Growth Medium
EC290	HiEncap™ Soyabean Casein Digest Agar
EC011	HiEncap™ Soyabean Casein Digest Medium
EC1316	HiEncap™ Super Broth
EC1689	HiEncap™ Super Broth No. II
ECG002	HiEncap™ Super Growth Agar
EC1250	HiEncap™ Terrific Broth (HiEncap™ Tartoff -Hobbs Broth)
ECG039	HiEncap™ YEP Growth Medium
EC1363	HiEncap™ YPD Broth (HiEncap™ YEPD Broth)
ECG038	HiEncap™ YPD Growth Agar (HiEncap™ YEPD Growth Agar)
EC1369	HiEncap™ YT Agar
EC1251	HiEncap™ YT Broth (HiEncap™ 2 XT Broth)
ECG031	HiEncap™ YT Growth Medium

ECXXXM
ECXXXD
ECXXXCCL

capsules contains media for preparing 1000 ml medium
capsules contains media for preparing 500 ml medium
capsules contains media for preparing 250 ml medium

HIMEDIA

Granulated Media

ONLY SINGLE COMPETITOR IN GLOBAL MARKET



Granulated Media

- More than 150 granulated medias listed in catalogue
- Developed for usersaerosol conscious
- Safe for use
- Less dust formation



Granulated Media

- Combines the high throughput technology of granulation and production of dehydrated culture media
- Has similar quality attributes with several added benefits in physical parameters



HIMEDIA®


For life is precious

Granulated Media

- Offering granulated media since 3 years
- Superior granulation technology
- Reduces aerosol formation
- Discrete flow properties
- Enhanced solubility
- Maximum wettability

HIMEDIA[®]

For life is precious



Granulated Media

- Granular forms can be manufactured for all existing media types
- Available for all applications
Food, Pharma, Clinical, Water, Cosmetics,
Dairy, Milk, Environment,
Research, Laboratory testing etc.

HIMEDIA[®]

For life is precious

Formulas & Testing Specifications

- Pharma guidelines
USP/BP/EP/JP/IP for pharma applications
- Food & Water as per ISO, FDA
BAM guidelines, APHA guidelines
- Clinical as per CLSI guidelines
Clinical handbook of Microbiology
- Other Standard Specifications

HIMEDIA[®]
For life is precious



Quality Control Testing*

Testing parameters

- Appearance
- pH
- Solubility
- Gelling
- Recovery (GPT Test)*
- Stability


HIMEDIA[®]

For life is precious



Diligent Quality Control

- Preparation of media
- Pouring of media
- Quantitative, semiquantitative and qualitative evaluations as in this standard
- Using Standard technique Results compared with standard reference medium (mandatory for quantitative methods) and compared with gold standard (previously approved batch)



HIMEDIA[®]

For life is precious

ISO 11133:2014

- Productivity & Selectivity limits
- Target & Non target organism
- Initial inoculum level defined for carrying out recovery testing
- PR of 0.7 ,PR of 0.5
- Recovery is $0.7 \times 100 = 70\%$, $0.5 \times 100 = 50\%$
- Quantitative testing and reporting on COA as specified in ISO 11133: 2014



For life is precious

HiMedia Laboratories Private Limited
23, Vadhani Industrial Estate, L.B.S. Marg, Mumbai - 400086
Website : www.himedialabs.com, Email : info@himedialabs.com

Certificate of Analysis, Quality and Conformity

Material Code : GM043	Material Name : Balrd Parker Agar Base Granulated	Lot No : 0000280394
Report No.: 040000683933	Date of Report : 11.11.2016	Expiry Date : Oct-2021

Appearance

Cream to yellow coloured, granular medium. Observed : Light yellow

Gelling

Firm, comparable with 2.0% agar gel.

Colour and Clarity of prepared medium

Basal medium: Yellow coloured clear to slightly opalescent gel. After addition of Egg Yolk Emulsion and Tellurite Emulsion: Yellow colour red opaque gel forms in Petri plates.

Reaction

Reaction of 6.3% w/v aqueous solution at 25°C.

pH

pH Range : 6.80-7.20 Observed : 7.10

Cultural response

Cultural response was observed after an incubation at 35-37°C for 24-48 hours. Recovery rate is considered as 100% for bacteria growth on Soyabe an Casein Digest Agar.

Organism	Inoculum (CFU)	Growth	Lot value (CFU)	Recovery	Colour of colony	Lecithinase
Cultural response						
<i>Staphylococcus aureus</i> ATCC 6538	87	luxuriant	82	94%	grey-black shiny	Positive, opaque zone around the colony
<i>Staphylococcus aureus</i> ATCC 25923	83	luxuriant	76	91%	grey-black shiny	Positive, opaque zone around the colony
<i>Proteus mirabilis</i> ATCC 25933	82	good - luxuriant	81	99%	brown - black	Negative
<i>Micrococcus luteus</i> ATCC 10240	80	poor - good	28	35%	shades of brown-black (very small)	Negative
<i>Staphylococcus epidermidis</i> ATCC 12228	89	poor - good	30	34%	black	Negative
<i>Bacillus subtilis</i> ATCC 6633	84	none - poor	3	3%	dark brown matt	Negative
<i>Escherichia coli</i> ATCC 8739	81	none- poor	6	7%	large brown black	Negative
<i>Escherichia coli</i> ATCC 25922	85	none- poor	6	7%	large brown black	Negative
<i>Escherichia coli</i> NCTC 9002	89	none- poor	4	4%	large brown black	Negative

- ATCC is a registered trade mark of the American Type Culture Collection

- NCTC and National Collection of Type Culture are registered trade mark of the Health Protection Agency

Control Media :

HiMedia Laboratories Private Limited

23, Vadhani Industrial Estate, L.B.S. Marg, Mumbai - 400086

Website : www.himedialabs.com, Email : info@himedialabs.com**Certificate of Analysis, Quality and Conformity**

Material Code : GM043	Material Name : Balrd Parker Agar Base Granulated	Lot No : 0000280394
Report No.: 040000683933	Date of Report : 11.11.2016	Expiry Date : Oct-2021

- For Bacteria : Soyabean Casein Digest Agar / Columbia Blood Agar base enriched with 5% v/v Sheep/Horse blood.
- For Yeast & Mold : Sabouraud Dextrose Agar.




- All ISO 11133 : 2014(E) control strains are included in the Quality parameter
- HiMedia Laboratories Pvt Ltd is certified for ISO 9001-2008, ISO 13485-2003 and WHO GMP.

- Information for BSE/TSE Risk: The material was subjected to pH \leq 7.0 and/or a temperature in excess of 75°C for no less than 2 hours during the manufacturing process. The bovine raw material for this product was collected entirely from Indian Origin animals in a licensed based establishment. The animals are inspected under a Govt. approved veterinarian's supervision and were apparently free from infectious and contagious diseases. BSE (Bovine Spongiform Encephalopathy)/ TSE (Transmissible Spongiform Encephalopathy) and dioxine are not known to exist in India. This material does not contain, nor is derived from the specific risks material as defined in The Maharashtra Animal Preservation Act Govt. of Maharashtra, India.

STATUS OF THE MATERIAL : APPROVED

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications . The information given here is believed to be correct and accurate, however, both the information and products are offered without warranty for any particular use, other than that specified in the current HiMedia manual or product sheets. The results reported were obtained at the time of release.

This document has been produced electronically and is valid

**Microbiologist/Analyst****Dy QC/Dy QA Manager****Quality Assurance Manager****11.11.2016**

Testing as per ISO 11133:2014

Cultural response

Cultural characteristics was observed after an incubation for Bacterial at 30-35°C 18-24 hours and for Fungal at 30-35°C <=5days.

Organism	Inoculum (CFU)	Observed Lot value (CFU)	Recovery	Observed Lot value (CFU) w/blood	Recovery w/ blood	Haemolysis
<i>Bacillus subtilis</i> ATCC 6633 (00003)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	none
<i>Staphylococcus aureus</i> ATCC 25923 (00034)*	50 -100	35 -100	>=70 %	35 -100	>=70%	beta
<i>Staphylococcus aureus</i> ATCC 6538 (00032)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	beta
<i>Escherichia coli</i> ATCC 25922 (00013)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	none
<i>Escherichia coli</i> ATCC 8739 (00012)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	none
<i>Escherichia coli</i> ATCC 11775 (00090)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	none
<i>Escherichia coli</i> NCTC 13167 (00179)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	none
<i>Escherichia coli</i> NCTC 9002	50 -100	35 -100	>=70 %	35 -100	>=70 %	none
<i>Pseudomonas aeruginosa</i> ATCC 27853 (00025)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	-
<i>Pseudomonas aeruginosa</i> ATCC 9027 (00026)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	-

HIMEDIA®

For life is precious

Testing as per ISO 11133:2014

<i>Pseudomonas aeruginosa</i> ATCC 10145 (00024)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	-
<i>Salmonella Abony</i> NCTC 6017 (00029)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	-
<i>Micrococcus luteus</i> ATCC 9341	50 -100	35 -100	>=70 %	35 -100	>=70 %	-
<i>Streptococcus pneumoniae</i> ATCC 6305	50 -100	35 -100	>=70 %	35 -100	>=70 %	-
<i>Salmonella Typhimurium</i> ATCC 14028 (00031)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	-
<i>Enterococcus faecalis</i> ATCC 29212 (00087)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	-
<i>Candida albicans</i> ATCC 10231 (00054)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	-
<i>Candida albicans</i> ATCC 2091 (00055)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	-
* <i>Aspergillus brasiliensis</i> ATCC 16404 (00053)*	50 -100	25 -70	50-70%			-
<i>Clostridium perfringenes</i> ATCC 13124 (00007)*	50 -100	35 -100	>=70 %	35 -100	>=70 %	-

Key : * - Corresponding WDCM numbers

HIMEDIA®

For life is precious

Labelling

- Meets Labelling as per harmonized standards EN980:2008 & ISO 15223-1:2012 for IVD and CE marking of products wherever applicable.
- Labelling as per CLP 1272/2008 for hazardous classification of products



For life is precious

Typical Label, Hazard



Storage :
Hygroscopic. Keep tightly closed
away from bright light. On receipt
store at



No liability accepted for accidents in handling or use

FOR LABORATORY USE ONLY

For In Vitro Diagnostics



HIMEDIA REF GM010-500G

Net Content 500g

Alternative Thioglycollate Medium, Granulated (Thioglycollate Broth Alternative, Granulated)

Directions:

Suspend 29.0 grams in 1000 ml purified / distilled water.
Heat if necessary to dissolve the medium completely.
Distribute into tubes or flasks as desired. Sterilize by
autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Note :

*It is preferable to use freshly prepared medium, alternatively
it should be boiled and cooled just once prior to use as on
reheating, toxic oxygen radicals are formed.*

*For More Information Refer Technical Data



Use

For sterility testing of turbid or viscous biological products as
per various pharmacopoeia and cultivation of anaerobes
from clinical and non-clinical specimens.

Warning

May cause an allergic skin reaction

Avoid breathing dust/fume/gas/mist/vapours/spray. Wear
protective gloves/protective clothing/eye protection/face
protection. IF ON SKIN: wash with plenty of soap and water.
IF SKIN irritation or rash occurs: Get medical
advice/attention.

**Standard Formula

Ingredients	Gms/litre
Tryptone	15.00
Yeast extract	5.00
Dextrose (Glucose)	5.50
Sodium chloride	2.50
L-Cystine	0.50
Sodium thioglycollate	0.50
Final pH (at 25°C)	7.1±0.2
**Formula adjusted, standardized to suit performance parameters	

COUNTRY OF ORIGIN-INDIA
Company certified for ISO 9001:2008,
ISO 13485:2003, WHO GMP



HiMedia Laboratories Pvt. Ltd.

Reg.off. : 23, Vadhani Ind. Est.,
LBS Marg, Mumbai-400086, India.

Works : B/4-6, M.I.D.C., Dindori, Nashik, India.

Customer Care No.: 00-91-22-6116 9797

Email: techhelp@himedialabs.com



MAR-2019



L0T 0000226786



8 902729 919330



Expected performance during specified expiry period when
material is duly maintained in the original powder form.

Received _____ Opened _____



*Disposal : User must ensure safe disposal by
autoclaving and/or by incineration used or unusable
preparations of this product and derivatives thereof
on completion of work to avoid contagion.*

HIMEDIA®

Typical Label, NonHazard



Storage :

Hygroscopic; Keep tightly closed
away from bright light. On receipt
store at



No liability accepted for accidents in handling or use

FOR LABORATORY USE ONLY

HIMEDIA[®] **REF** **GMH290-500G**

Net Content 500g

☉ Equivalent to MU290/ME290/MM290/M290B

Soybean-Casein Digest Agar (Casein Soyabean Digest Agar), Granulated

Directions:

Suspend 40.0 grams in 1000 ml purified/
distilled water. Heat to boiling to dissolve the
medium completely. Sterilize by autoclaving
at 15lbs pressure (121°C) for 15 minutes or
as per validated cycle. Cool to 45-50°C. Mix
well and pour into sterile Petri plates.

*For More Information Refer Technical Data

Use

Recommended as a general purpose medium for cultivation
of a wide variety of microorganisms from pharmaceutical
products in accordance with the harmonized method of
USP/EP/BP/JP/IP (Medium 2).

**Standard Formula

Ingredients	Gms/litre
Tryptone \$	15.00
Soya peptone [^]	5.00
Sodium chloride	5.00
Agar	15.00

*pH after sterilization 7.3±0.2

*pH can also be measured after sterilization at 25°C

**Formula adjusted, standardized to suit performance
parameters

\$ Pancreatic digest of casein

[^]Papaic digest of soyabean meal

COUNTRY OF ORIGIN-INDIA

Company certified for ISO 9001:2008,
ISO 13485:2003, WHO GMP

HiMedia Laboratories Pvt. Ltd.

Reg.off. : 23, Vadhani Ind. Est.,
LBS Marg, Mumbai-400086, India.

Works : B/4-6, M.I.D.C., Dindori, Nashik, India.

Customer Care No.: 00-91-22-6116 9797

Email: techhelp@himedialabs.com



FEB-2022



LOT 0000292057



8 902729 914632



Received _____ Opened _____



Disposal : User must ensure safe disposal by
autoclaving and/or by incineration used or unusable
preparations of this product and derivatives thereof
on completion of work to avoid contagion.

Expected performance during specified expiry period when
material is duly maintained in the original powder form.

HIMEDIA[®]

For life is precious

Old Label, Hazard



Granulated

HIMEDIA



GM009-500G

Net Content
500g

Fluid Thioglycollate Medium, Granulated (Thioglycollate Medium, Fluid, Granulated)

It is used for stability testing of biologicals and for cultivation of anaerobic, aerobic and microaerophilic

FOR LABORATORY USE ONLY

Storage :

Hygroscopic. Keep tightly closed, away from bright light. On receipt store in a cool, dry place. It is not to be accepted for use unless it is found to be free from any contamination.

COUNTRY OF ORIGIN-INDIA

Directions :

Suspend 29.75 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 25°C and store in a

cool dark place preferably below 25°C.

Note : If more than the upper limit of the medium has evaporated during storage, the medium may be restored once by heating in a water bath or in the flowing steam until the pink colour disappears.



NO
GHS

NO
GHS

Warning

May cause an allergic skin reaction. Avoid breathing dust/fume/gas/mist/vapour/spray. Wear protective gloves/protective clothing/eye protection/face protection. FOR SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.

*Standard Formula

Ingredients	Gm/Litre
Tryptone	15.00
Yeast extract	5.00
Dextrose (Glucose)	5.50
Sodium chloride	2.50
L-Cystine	0.50
Sodium thioglycollate	
Final pH (at 25°C) 7.1 ± 0.2	
*Formula adjusted standardized to suit performance parameters.	

Company certified for ISO 9001:2008, ISO 13485:2003, WHO GMP

Himedia Laboratories Pvt. Ltd.

Works : B-4-B, M.I.D.C., Dindori, Nashik, MH.
Regd. Off. : 23, Vaidya Hall Ind. Est.,
LBS Marg, Mumbai-400006, India.
Customer Care No. : 00-91-22-61116 9797
Email : techhelp@himedialabs.com



EXP-2018



LOT 00000020627



9 782729 114258



Received: _____ Opened: _____



Disposal : Use most appropriate disposal by autoclaving and/or by incineration and/or suitable preparation of this product and derivatives thereof on completion of work to avoid contamination.

Inspect performance during specified expiry period when water test is performed in the original container.

EXP-2018

MUM/PC/Pha/200004- Incl. of all labels

HIMEDIA

For life is precious

Safety data sheet



www.himedialabs.com

Safety data sheet(SDS)

According to Regulation (EC) No.1907/2006

Revision : 00001

Date of Revision : 09.01.2017

1 Identification of the substances/ mixture and of the company/ undertaking

1.1 Product Identifiers

Product Number GM024
Product Name Cetrimide Agar Base, Granulated
REACH Registration Number This product is a mixture. Reach registration number is not available for this mixture.

1.2 Relevant identified uses of the substance or mixture and uses advised against

1.2.1 Relevant identified uses Laboratory Chemicals, Analytical Purpose, Biochemical Analysis
For InVitro Diagnostic Use

1.3 Details of the supplier of the safety data sheet

Produced by HiMedia Laboratories Private Limited
Address 23, Vadhani Industrial Estate, LBS Marg, Ghatkopar (W), Mumbai - 400 086
India
Tel. No. +91-22-2500 0970, +91-22-2500 1607 Fax No. : +91-22-25002468
Mail Id info@himedialabs.com Website : www.himedialabs.com

1.4 Emergency Tel. No.

Emergency Tel. No. Please contact the regional HiMedia representation in your country

Safety data sheet

Component		Classification	Concentration
Cetrimide			
CAS No. :	57-09-0	As Per EC Regulation 1272/2008 Acute Tox.oral 4; Skin Irrit. 2; Eye Irrit. 2A; STOT SE 3; Aquatic Chronic 1 H302; H315; H319; H335; H410	>=0.1 - <=1.0%
EC No. :	200-311-3		

Refer Section 16 for complete statement of H codes & classification.

4 First Aid Measures

4.1 Description of first aid measures

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

In case of skin contact

Wash with plenty of soap and water. Consult a physician.

In case of eye contact

Rinse immediately with plenty of water for at least 15 minutes. Consult a physician.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

4.2 Most important symptoms and effects, both acute and delayed

No data available.

4.3 Indication of immediate medical attention and special treatment needed

No data available

Safety data sheet

2 Hazards Identification

2.1 Classification of the substance or mixture

CLP Classification-Regulation (EC) No. 1272/2008[EU-GHS/CLP]

Hazardous to the aquatic environment, long term hazard, (Category 1), H410

2.2 Label elements

Labeling according to Regulation (EC) No.1272/2008



Pictogram

Signal word Warning

Hazard Statement(s)

H410 Very toxic to aquatic life with long lasting effects

Precautionary Statement(s)

P273 Avoid release to the environment.

2.3 Other Hazards

None

3 Composition/Information On Ingredients

3.2 Mixture

Safety data sheet

5 Fire Fighting Measures

5.1 Extinguishing media

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Unsuitable extinguishing media

No data available.

5.2 Special hazards arising from the substance or mixture

Carbon oxides, Sulphur oxides, Hydrogen chloride gas, Potassium oxides, Magnesium oxide

5.3 Precautions for fire-fighters

Wear self contained breathing apparatus for fire fighting if necessary

5.4 Further information

No data available

6 Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures



THANK YOU

Any Questions Please ???
Or write to
hdmexpo@himedialabs.com
for further queries



For life is precious